

LU-CB1 Leiden

Micro-Incubator

User's Manual

LU-CB1 Leiden Micro-Incubator User's Manual 65-0040



HARVARD
A P P A R A T U S

Table of Contents

1

LU-CB1 Leiden Micro-Incubator User's Manual

SUBJECT	PAGE
Warranty and Repair Information	2
Serial Numbers	2
Calibrations	2
Warranty	2
Repair Facilities and Parts	2
Introduction	3
Micro-Incubator Construction	3
Operation	4
Perfusion	5
Fluid Perfusion	5
Gas Perfusion and Grounding	6
Gas Perfusion	6
Grounding	6
Advantages and Accessories	7
Included Accessories	7
Optional Accessories.....	7
Setup	8

Warranty and Repair Information

Serial Numbers

All inquiries concerning our product should refer to the serial number of the unit. Serial numbers are located on the rear of the chassis.

Calibrations

All electrical apparatus is calibrated at rated voltage and frequency.

Warranty

Harvard Apparatus warranties this instrument for a period of one year from date of purchase. At its option, Harvard Apparatus will repair or replace the unit if it is found to be defective as to workmanship or material.

This warranty does not extend to damage resulting from misuse, neglect or abuse, normal wear and tear, or accident.

This warranty extends only to the original customer purchaser.

IN NO EVENT SHALL HARVARD APPARATUS BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow exclusion or limitation of incidental or consequential damages so the above limitation or exclusion may not apply to you. **THERE ARE NO IMPLIED WARRANTIES OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR USE, OR OF ANY OTHER NATURE.** Some states do not allow this limitation on an implied warranty, so the above limitation may not apply to you.

If a defect arises within the one-year warranty period, promptly contact ***Harvard Apparatus, Inc. 84 October Hill Road, Building 7, Holliston, Massachusetts 01746-1371*** using our toll free number 1-800-272-2775. Goods will not be accepted for return unless an RMA (returned materials authorization) number has been issued by our customer service department. The customer is responsible for shipping charges. Please allow a reasonable period of time for completion of repairs, replacement and return. If the unit is replaced, the replacement unit is covered only for the remainder of the original warranty period dating from the purchase of the original device.

This warranty gives you specific rights, and you may also have other rights which vary from state to state.

Repair Facilities and Parts

Harvard Apparatus stocks replacement and repair parts. When ordering, please describe parts as completely as possible, preferably using our part numbers. If practical, enclose a sample or drawing. We offer a complete reconditioning service.

CAUTION


This apparatus is not registered with the FDA and is not for clinical use on human patients.



CAUTION: Not for clinical use on human patients.

Introduction

3



The LU-CB1 is a simple micro-CO₂-incubator designed for use on the microscope stage of an inverted microscope. The micro-incubator is easy to use and offers a handy tool for the culture of cells under the microscope. The performance of the LU-CB1 compares well with that of a conventional CO₂-microincubator.


A standard 35 mm plastic or glass culture dish can be placed in the micro-incubator. The culture medium is then covered by a gas permeable layer of mineral oil, which protects the culture from the environment without affecting culture conditions, allowing easy cell manipulation under microscopical control.

Micro-Incubator Construction

The LU-CB1 Leiden Micro-Incubator consists of:

- 3 torroidal rings anodized with black aluminum-oxide film
- Two-turns of cantal or nichrome wire wound around the middle ring and connected to external 1 mm jacks
- Gas perfusion inlet
- Polyvinylchloride (PVC) support ring with magnetic base (metal top) to allow mounting of thermistor probes or other devices with magnetic holders around the incubator. (Support ring is fitted with wells that can be used as reference electrode wells for electrophysiological experiments and two screws to hold the micro-incubator in place.)
- 3 (PVC) insulating rings to avoid contact with the insert plate of the microscope stage and to prevent heat loss.

Operation



The incubator is heated by passing electric current through the heating coil inside the anodized aluminium ring. The duct inlet in the ring directs and heats the incoming gas. The cell culture is thus heated by direct heat conduction from the incubator as well as by convection as result of the heated gas blown over the surface of the bath or chamber. The heat transfer by convection reduces considerably the temperature gradient between the center and side of the culture dish. A 35 mm OD x 10 mm high disposable culture dish or a Glass cover slip dish fits snugly into the incubator. A layer of non-toxic mineral oil (2 to 3 ml) added to the dish will prevent the evaporation of the medium in the dish (2 to 4 ml).

For temperature control, the TC-202A and its appropriate sensor BSC T3 are used. The temperature control unit regulates the heating current by sensing the temperature of the incubator or of the medium via its sensor.

The TC-202A has been designed as bi-polar and mono-polar unit. To achieve the feed-back control and obtain maximum benefits from your LU CB1 Micro-Incubator, the TC-202A must be used on its mono-polar mode. (See TC-202A setup for more information.)

The TC-202A has been designed as bi-polar and mono-polar unit. To achieve the feed-back control and obtain maximum benefits from your LU CB1 Micro-Incubator, the TC-202A must be used on its mono-polar mode. (See TC-202A setup for more information.)





Figure 1. Cross-Sectional Drawing of Micro-Incubator




Fluid Perfusion

The perfusion system above uses a peristaltic pump for both inflow and outflow of the bath. On entry into the Faraday cage, a bubble trap ensures electrical isolation of each channel. The volume of fluid in the chamber is regulated by the changing height of the suction tube. (LU ASP) Although the micro-incubator was not designed with internal perfusion ports, perfusion can be achieved by using a system that allows the fluid to be pre-heated before entering in the bath. However, oil cannot be used in this configuration.

The temperature of the medium and the pH should be controlled just before it enters the culture dish. The fluid temperature must also be as close as possible to the temperature in the culture dish and have a constant flow rate, minimizing the temperature oscillation of the medium. The heating system of the culture bath then can be used to support the pH and the temperature control of the medium as it passes through the culture dish.

For reproducibility at low flow rates, a peristaltic pump is preferable over gravity feed. This configuration also makes it easier to maintain the fluid level while changing solutions. An aspirator (LU ASP) can be used to remove the perfusate in the chamber via direct suction. The aspirator has a magnetic base to grip the magnetic metal base on top of the support ring. The level of fluid in the chamber is determined by the height of this aspirator and can be changed by adjusting the nut on top of the threaded post. The oscillations inherent in any peristaltic perfusion are easily damped with bubble-traps.



Gas Perfusion and Grounding

6

Gas Perfusion

This technique is intended primarily to reduce the radiant heat loss to (or gain from) the local atmosphere, thereby reducing temperature gradients across the dish by controlling the gas mixture above the chamber. The gas passes over the heat exchange. Any desirable gas mixture may be introduced into the micro-incubator.

Besides its use to achieve the lowest possible temperature gradients between the center as the side of the culture dish, the gas can be used to maintain the partial pressure of gases in the culture medium (e.g., 5% CO₂ and 95% at a flow rate of 0.5 to 2 L/min, the pH of a bicarbonate-buffered medium and its pO₂ can be maintained).

Grounding

Two independent ground pathways are also available when single-channel or similar demanding electrical recordings are to be performed.

1. **Salt Bridge Ground:** The Support ring is fitted with a well that can be used as Bath Earth Pot for electrophysiology experiments. A silver/silver chloride (Ag/AgCl) electrode can be introduced in the pot, then is filled with KCl solution and an AGAR bridge used to complete the circuit with the medium in the culture dish during the experiment.
2. **Shield Ground:** A grounding terminal is connected at the end of the output cable, connecting this to recording instruments or local grounds such as microscopes or Faraday cage.

Advantages and Accessories

7

The two main advantages of the oil-covered micro-incubator are that it is freely accessible for experimental manipulations and that disposable or reusable culture dishes can be used. For biochemical analysis, aliquots can be taken from the culture under the mineral oil surface, and medium can be easily refreshed or exchanged in the same way.


In light microscopical studies, for example, a transparent lid can be placed on the micro-incubator to reduce the gas flow and drugs can be easily applied to the culture by opening and closing this lid. During electro-physiological experiments the lid can not be used, however, because free access to microelectrodes is necessary. Gas flow then should be at least 1 liter/min. to minimize temperature gradient and obtain adequate equilibration of the gas and the medium. For all cases the use of more than 2 ml medium is recommended.

Included Accessories

1. Transparent Cover Lid
2. Pin Magnet
3. Set of 3 Insulating Rings with 40 mm, 48 mm and 50 mm Diameters
4. 6 ft. Output Cable
5. Sensor or Thermistor Holder
6. 1 mm Pins Spare (2 pcs.)

Optional Accessories

1. Peristaltic Pump (MP II)
2. Aspirator (LU ASP) for Fluid Perfusion
3. Alignment Rings to replace Insert Microscope Stage for Nikon, Olympus, Zeiss or Leica Microscopes
4. Fixing Platform (LU FPZ) for Zeiss Axiovert with Attachable Mechanical Stage, M Type 451335/6
5. Fixing Platform (LU_FPL) for Leica DMIRB Mechanical Attachable Stage 520804
6. Cover Slip Glass Teflon[®] Dish



1. Set the temperature controller to Mono-polar. (This is done by setting the Incubator selector switch in the back of the TC-202A to: "BSC-BU; LU-CB1 & LU-CPC".)
 2. Assemble the micro-incubator into the appropriate support ring as follows:
 - a. Align the 1 mm sockets and the CO₂ duct with the cutout of the support ring.
 - b. Insert the micro-incubator all the way up, starting at the bottom of the support ring.
 - c. Secure the micro-incubator by tightening the screws located in the opposite side of the cutout of the support ring.
 - d. Assemble one of the PVC insulator rings (depending on the size of the opening of the stage insert plate) to the outside of bottom of the micro-incubator. This ring will act as an insulator to prevent heat loss by insulating the LU-CB1 from the microscope stage.
3. Secure the bottom section of the support ring to the stage insert plate or PDMI Alignment Ring (by using a sticky compound such as modeling clay, Silicon Glue or similar method).
 4. Connect the 6 pin Din connector output cable (Catalog No. 65-0100) to the TC-202A Front Panel (I/O).
 5. Insert the cables with the 1mm plug through the strength relief holes and connect the 1mm pin to the sockets of the micro-incubator. The ring terminal can be used to ground any other instrument related to your experiment.
 6. Connect the temperature probe to the Bath Sensor in the TC-202A.
 7. Insert the 35 mm culture dish into the micro-incubator.
 8. Attach the magnetic thermistor holder on the magnetic base (metal plate on top of the support ring). Insert the temperature/sensor probe through the coil and into the bath. The temperature will be displayed in the LCD of the TC-202A.
 9. Set the "Control Point" of the TC-202A to "BATH."
 10. Turn the TC-202A "ON" and adjust the set point to the desired temperature (from ambient up to 50⁰C).
 11. Depending on the ambient the temperature, the bath temperature will stabilize within 15 to 20 min after initial set up.