


# microelectrode holders

## Electrode Holder Replacement Parts


### Microelectrode Replacement Parts

Order #	Qty/Pkg	Model	Description
<b>Coupling Elements</b>			
W4 64-1297	3	SP-1M	Molded Ag-AgCl pellet assembly
W4 64-1282	10	AG25-10	Ag wire, 0.25 mm diameter x 70 mm long
W4 64-1319		AG10W	Ag wire, 0.25 mm diameter x 2 m (6'6")
W4 64-1288	1	HP-205	Hybrid Ag-AgCl pellet/Ag wire assembly with 40 mm teflon sleeve sealed at junction of pellet and silver wire




### Pipette and Wire Seals

W4 64-1289	10	PS-10	Pipette seal for 1.0 mm diameter glass
W4 64-1290	10	PS-12	Pipette seal for 1.2 mm diameter glass
W4 64-1291	10	PS-15	Pipette seal for 1.5 mm diameter glass
W4 64-1374	10	PS-17	Pipette seal for 1.7 mm diameter glass
W4 64-1292	10	PS-20	Pipette seal for 2.0 mm diameter glass
W4 64-1298	10	WS-1	Wire seal for E, ME, MP and theta holder
W4 64-1299	10	WS-2	Wire seal for Q and PE Series holders
W4 64-1300	10	WS-3	Wire seal for Q and PE Series holders with threaded "T" connector



### Glass Seal Compression Caps

W4 64-1293	2	QC-10	Compression cap for 1.0 mm glass holders
W4 64-1294	2	QC-12	Compression cap for 1.2 mm glass holders
W4 64-1295	2	QC-15	Compression cap for 1.5 mm glass holders
W4 64-1375	2	QC-17	Compression cap for 1.7 mm glass holders
W4 64-1296	2	QC-20	Compression cap for 2.0 mm glass holders



### Connector Pins and Jacks

W4 64-1283	3	HC-10M	1 mm pin for Q and PE holders with A or T connectors
W4 64-1284	3	HC-13M	1.3 mm pin for Q and PE holders with B (BNC) connector
W4 64-1285	3	HC-20M	2 mm threaded pin for E, ME and theta holders
W4 64-1286	3	HC-21F	2 mm jack assembly for all series holders
W4 64-1287	3	HC-22M	2 mm pin for Q and PE holders with M connector