



Plugs and Jacks

An assortment of gold-plated 1 mm plugs and jacks, with or without wires. Useful for grounding, shielding, etc.

Order #	Model	Product
W4 64-1325	WC1-10	1 mm Pin, pkg. of 10
W4 64-1326	WC2-5	2 mm Pin, pkg. of 5
W4 64-1322	PJ1-10	1 mm Jack, Uninsulated, pkg. of 10
W4 64-1323	PJ2-5	2 mm Jack, Uninsulated, pkg. of 5
W4 64-1324	TJ1-3	1 mm Jack, Teflon® Insulated, pkg. of 3
W4 64-1328	WA30-5	1 mm Pin with 30 cm, 26 ga. Insulated Wire, pkg. of 3
W4 64-1327	WA10-5	1 mm Pin with 10 cm, .25 mm D, Bare Silver Wire, pkg. of 2 For Bilayer Headstage

Silver Wire

Made of annealed silver 99.99% pure. Supplied for the convenience of researchers to fashion their own electrodes.

Order #	Model	Product
W4 64-1318	AG8W	Silver Wire, 2 m, 0.20 mm (.008") D, L, No Insulation
W4 64-1319	AG10W	Silver Wire, 2 m, 0.25 mm (.010") D, L, No Insulation
W4 64-1320	AG15W	Silver Wire, 2 m, 0.37 mm (.015") D, L, No Insulation
W4 64-1321	AGT10W	Silver Wire, 2 m, 0.25 mm (.010") D, L, Teflon® Insulation