

# bilayer workstation

## Additional Components

RAC-14



### RAC-14 Instrument Rack

An instrument rack provides an organizing influence on the electronics and facilitates the tabletop design of the Workstation.

The RAC-14 is a 14" tabletop rack and comes with all the hardware needed to mount your amplifier, filter, stirring-illumination controller and A/D board into one integrated unit. A ground point is provided for attachment to a Power Line Conditioner. Powder black finish.

Order #	Model	Product
W4 64-0070	RAC-14	Instrument Rack

### BLM Starter Kit

Warner Instruments has put together a Starter Kit to aid you in completing your Bilayer Workstation. The kit includes 3 BNC cables, grounding cables, a Red Sable artist's dotting brush, glass capillary tubing, PE tubing (for your perfusion system), and one replacement pack of stirbars (2 x 5 mm).

Order #	Model	Product
W4 64-0067	BLM-ST	Bilayer Starter Kit

HST-1 Headstage Holder Kit



### HST-1 Headstage Holder Kit

- Shock mounting for headstage
- Non-conducting surface
- Flexible headstage attachment
- Magnetic base with on/off switch

The HST-1 is designed as a shock mount positioning device for the bilayer headstage. The non-conducting neoprene pad, mounted onto a polycarbonate base, serves as a stable platform for the headstage. Two medium elasticity rubber bands are used to fix the headstage to the pad and provides a secure but flexible attachment.

Convenient positioning of the HST-1 Headstage Holder is provided by an MB/B magnetic base. The base contains a magnet with on/off switch for easy re-positioning of the assembly within your application.

### Specifications

Platform, H x W x D	5.7 x 4.4 x 1.3 cm (2.24 x 1.12 x 0.5 in)
Material	Polycarbonate base; non-conducting neoprene pad
Shaft, Dia x L	0.6 x 6.4 cm; anodized aluminum

Order #	Model	Product
W4 64-0435	HST-1	Headstage Holder System

### Replacement parts

W4 64-0422	HH-1	Headstage Holder
W4 64-0060	MB/B	Magnetic Base