

culture dish platforms

Accessories

Chambers

culture dish platforms

Heated platforms provide a base for Series 20 chambers and allow easy clamping of coverslip and chamber

Series 20 Interface Cables and Thermistors



CC-15



CC-28



CC-35



TA-29




TA-30



TS-60P

1.95 mm x 9 mm



TS-70B

0.85 mm x 3.6 mm



RH-2

Order #	Model	Product
W4 64-0303	CC-15	Interface Cable with connector with tinned leads
W4 64-0106	CC-28	Interface Cable with connector for TC-324B/344B
W4 64-0109	CC-35	Cable with with 8 pin DIN connector for TC-324B/344B
W4 64-0107	TA-29	Cable with Bead Thermistor

Order #	Model	Product
W4 64-0108	TA-30	Cable with Glass Bead Thermistor
W4 64-0269	TS-60P	Replacement Probe Thermistor for CC-28 and CC-15 Cables
W4 64-0270	TS-70B	Replacement Bead Thermistor for CC-28 and TA-29 Cables
W4 64-0274	RH-2	Replacement Heater Elements for PH and PM Platforms

Platform and Chamber Parts

Order #	Model	Product
W4 64-1429	STU-1	Suction Tube Upgrade Kit with Mag-7
W4 64-1506	RPC-1	Replacement Platform Clamps
W4 64-1507	RPC-2	Replacement Clamps for DH35i, DH-40i and QE-1
W4 64-0276	CR-15	Coverslip Retainer for RC-20 and RC-20H



STU-1



RPC-1 Clamps & Screws



RPC-2



CR-15

culture dish platforms

Accessories (continued)

MCK-1 Magnetic Clamp Kit for Series 20 Chambers




- Useful for simple positioning needs
- Complete kits available for all Series 20 platforms
- Includes a variety of holding tools

The MCK-1 Magnetic Clamp Kit provides our classic Series 20 chamber users with a convenient method of mounting or clamping accessories to the perfusion chamber. Typical items that may be attached are additional suction or perfusion tubes, ground electrodes, hold down needles, thermistor probes or gas jets.

Each kit comes with 4 sets of magnetic stainless steel mounting plates that fit on all P and PH platforms and 6 magnetic supports. A clamp for mounting 1/4 inch (6.3 mm) diameter devices has a rotating joint with locking knob. This clamp is useful for turning a standard electrode holder into an agar bridge for solution bath grounding. Two included low profile twin magnets with silicone rubber pads can be used for mounting needle holddowns. The kit also includes a mechanical clamp with two slits that can clamp onto small diameter tubing such as a suction tube. The remaining two clamps are designed for securing ground wires or thermistors.

Order #	Model	Product
W4 64-0357	MCK-1	Magnetic Clamp Kit
W4 64-0358	MAG-1	Magnetic Clamp for 1/4 inch (6.3 mm) devices
W4 64-0359	MAG-2	Magnetic Clamp 16 g for suction tubes
W4 64-0360	MAG-3	Magnetic Clamp for wires
W4 64-0361	MAG-4	Magnetic Clamp with twin magnetics



Made from magnetic stainless steel and delrin this kit is easily added to your set-up



Each kit comes with four complete base plate designs, offering considerable versatility.



Magnets with just the right amount of force are used in the clamp kit



Suction tube holder