

Coverslips chamber accessories

Large selection of high quality coverslips stocked for immediate delivery



Now Available in #1.5 Thickness!

High quality Coverslips are essential for microscopy imaging. At high resolutions (40X, 60X or 100X), choosing the correct coverslip is critical for good performance. Warner has recently expanded its line of coverslips to include German borosilicate glass at #1.5 thickness. Except CS-40R, CS-30R and CS-24/60.

Stocked in depth in the sizes used with RC Series, Series 20 and the new Series 40 chambers. Consult the chart to select the correct cover glass for your chamber.

Order #	Model	Quantity	Dimension	Chamber Used With
Coverslip Thickness #1: 0.15 mm (0.006 in.)				
W4 64-0720	CS-3R	100	3 mm D	
W4 64-0700	CS-5R	100	5 mm D	Fits inside 96-well culture plate
W4 64-0701	CS-8R	100	8 mm D	Fits inside 24-well culture plate
W4 64-0702	CS-12R	100	12 mm D	RC-25, RC-48LP
W4 64-0703	CS-15R	100	15 mm D	RC-20, RC-20H, RC-25F, RC-37WC, RC-37FC
W4 64-0384	CS-18R	100	18 mm D	RC-37WC, RC-41LP, RC-43C
W4 64-0704	CS-22S	100	22 x 22 mm	RC-21B
W4 64-0705	CS-25R	100	25 mm D	RC-21BR, RC-21BRFS, RC-21BRW, RC-21BDW, RC-40LP, RC-40HP, RC-43C, RC-46SNLP, RC-46SLP, RC-47FSLP, RC-50
W4 64-1500	CS-40R	60	40 mm D	RC-30WA
W4 64-0706	CS-22/30	50	22 x 30 mm	RC-30, RC-30HV
W4 64-0707	CS-22/40	50	22 x 40 mm	RC-22, RC-22C, RC-24, RC-24E, RC-24N, RC-26, RC-26G, RC-26GLP, RC-26GS, RC-30, RC-50
W4 64-0708	CS-22/50	50	22 x 50 mm	RC-26Z, JG-23N/HP, JG-23W/HP, JG-23W/LP
W4 64-0709	CS-24/50	50	24 x 50 mm	RC-27, RC-27L, RC-27N, RC-27NE, RC-28, RC-29
Coverslip Thickness #1.5: 0.17 mm (0.0067 in.)				
W4 64-0718	CS-10R15	100	10 mm D	
W4 64-0712	CS-12R15	100	12 mm D	RC-25, RC-48LP
W4 64-0713	CS-15R15	100	15 mm D	RC-20, RC-20H, RC-25F, RC-42LP, RC-37WC
W4 64-0714	CS-18R15	100	18 mm D	RC-37FC, RC-41LP, RC-43C
W4 64-0719	CS-22R	100	22 mm D (Round)	
W4 64-0721	CS-22S15	100	22 mm D (Square)	
W4 64-0715	CS-25R15	100	25 mm D	RC-21BR, RC-21BRFS, RC-21BRW, RC-21BDW, RC-40LP, RC-40HP, RC-43C, RC-46SNLP, RC-46SLP, RC-47FSLP, RC-50
W4 64-1499	CS-30R	90	30 mm D	RC-30WA
W4 64-0716	CS-22/3015	50	22 x 30 mm	
W4 64-0717	CS-22/4015	50	22 x 40 mm	RC-22, RC-22C, RC-24, RC-24E, RC-24N, RC-26, RC-26G, RC-26GLP
Coverslip Thickness #2: 0.22 mm (0.0086 in.)				
W4 64-0710	CS-24/60	40	24 x 60 mm	RC-10, RC-11, RC-13, RC-16, RC-27D, RC-27LD