

## Temperature Control Accessories and Cables

### SHH-1, SHH-2, SHH-3 and SHH-4 Holders for In-line Solution Heaters



An ideal tool for holding Warner In-line Solution Heaters close to the chamber. Machined from a solid Delrin block, these sturdy holders can be secured to your microscope stage using gaffer or duct tape.

- The SHH-1 works with both our SH-27B and SF-28 Solution Heaters.
- The SHH-2 mounts our SC-20 In-line Solution Heater/Cooler.
- The SHH-3 is used for the SHM-6, SHM-8, and SHM-828 Multiline In-line Heaters.
- The SHH-4 works with our FR-50 and FR-55S flow valves.

Order #	Model	Product
W4 64-1555	SHH-1	Holder for Solution Heaters Models SH-27B & SF-28
W4 64-1556	SHH-2	Holder for Solution Heaters Model SC-20
W4 64-1557	SHH-3	Holder for Solution Heaters Models SHM-6, SHM-8, SHM-828
W4 64-1558	SHH-4	Holder for Flow Valves Models FR-50 and FR-55S

### Cables and Thermistors



Order #	Model	Product
W4 64-0303	CC-15	Cable Assembly with Connector with Tinned Leads
W4 64-0106	CC-28	Cable Assembly for Series 20 Heater Platforms
W4 64-0109	CC-35	Cable Assembly with Unterminated Outboard End
W4 64-1425	CC-102	Cable Assembly, Male Banana to 1 mm Jacks (for use with Warner Field Stimulation Chambers)
W4 64-0107	TA-29	Replacement Cable with Bead Thermistor
W4 64-0108	TA-30	Replacement Cable with Glass Thermistor

Additional replacement parts are shown on page 62 (Platforms).